

Conflicts of Interest

Consultant:
Microvention,
Medtronic, Stryker,
Boston Scientific

Scientific Board Member: Peytant







6 SIGMA OF ANEURYSM COILING

The 6 Sigma of Coiling

- 1. Size
- 2. Shape
- 3. Stiffness
- 4. Surface Modification
- 5. SAH
- 6. Support (Microcatheter/Balloon/Stent)

HOW TO CHOOSE THE FRAMING/FIRST COIL





- Choose the average diameter.
- 3D COIL.
- Choose the longest coil with that diameter that will fit, then decice on something slightly shorter.

- SIZE FOR THE LARGER DIMENSION.
- 3D COIL; BOXY SHAPE.
- GO RELATIVELY FAST WHEN PLACING ACROSS NECK
- PLAN TO USE A BALLOON

- SIZE FOR THE SMALLER DIAMETER.
- HELICAL IS A GOOD CHOICE.
- THE LONGEST COIL, WHICH CAN FIT BUT USE A SLIGHTLY SHORTER COIL.
- POSITION
 MICROCATHETER
 DEEP

- STENT.
- 3D COIL.
- LARGEST DIAMETER.
- DON'T GO TOO LONG.
- JAIL IF POSSIBLE.





- 3D COIL
- COMPLEX SHAPE SOFT
- DIAMETER ADJUSTED TO NECK AND BODY OF ANEURYSM.
- DO NOT GO TOO LONG

SIZE AND SHAPE: THE SECOND COIL

THERE IS NO RULE ON DOWNSIZING

- Need not always be helical
- Need not always be one size lower


THE DECISION
SHOULD BE
BASED ON
FINISH OF FIRST
COIL AND
APPEARANCE

DO NOT PUSH IT, PAINT IT.

KEEP WITHIN THE LANES


SIZE AND SHAPE: THE LAST COIL

- GO SHORT; FILL WITH PEBBLES NOT ROCKS
- THE ENEMY OF GOOD IS PERFECT
- COMPACTION WONT KILL, COMPLICATIONS WILL
- WHEN YOU WANT ONE MORE... DO NOT HAVE IT
- GEL IF YOU CAN
- GO SOFT
- DON'T GO SO SHORT THAT IT FALLS OUT



$$k = \frac{D_1^4 G}{8D_2^3 n} = Stiffness \propto \frac{D_1 G}{D_2 n}$$
.





Bending stiffness =
$$\frac{GD_1^4 p \cos(\gamma)}{16D_2 \left[\frac{1 + \sin^2(\gamma)}{1 + \nu} + \cos^2(\gamma) \right]}$$
$$= \frac{GD_1^4 \pi \sin(\gamma)}{16 \left[\frac{1 + \sin^2(\gamma)}{1 + \nu} + \cos^2(\gamma) \right]}$$




STIFFNESS AND SELECTION: FRAMING

STIFFNESS AND SELECTION: FILLING


STIFFNESS AND SELECTION: FINISHING

Surface Modification







Histology




- Softer coils are better for most aneurysms
- Balloon assistance can help reduce intra-operative rupture but can also have higher thromboembolic complications
- Packing density is important
- Give hydrocoils time to swell: interstices might stop filling
- Staged procedures are not uncommon






CASE SERIES: SAH



83 year old with HH3 SAH






Headway Duo 156, Scepter XC 4 x 11







Penumbra extra soft, hydrosoft, Penumbra wave






43 Year old with HH 4 SAH






Headway 17, Scepter XC 4 x 11; Hydrocoils





FBTS









LVIS JR: 2.5 x 23: keep calm and coil on

Final result







50 year old with hemiparesis



DSA




Balloon remodeling




- Scepter XC 4 x 11 mm
- Hydrocoils

Post embolization




Coiling through the Scepter

Coiling through the Scepter





54 y/o with SAH
















60 year old with SAH



Sofia and Scepter XC






Post coiling




66 year old with SAH






Scepter XC and Traxcess






Post coiling






45 y/o with SAH







Scepter XC 4 x 11



Stent coil with LVIS Blue ..45 year old with SAH



Extremely irregular PCOM aneurysm






Initial Coil Embolization: Scepter XC






Delayed finish: LVIS Blue 3.5 x 22, Duo Microcatheter for coiling.





Final Construct






Excellent wall apposition in a tight turn






Final DSA result





Simultaneous coiling: 69 year old with ruptured giant aneurysm at the top of the basilar $(2.9 \times 2.6 \times 2.6 \times 1.6 \times 2.6 \times 1.6 \times 1.6$





69 year old with ruptured giant top of basilar aneurysm 2.9 x 2.6 x 2.6 cm in size







It was difficult to visualize the PCAs on single vessel injection, a 7 french guide cath in the left vert and a 5 french guide cath in the right vert both had to be injected for 14 seconds before start of DSA to see PCAs, 3D runs also needed a 14 second delay







2 Scepter XC catheters were advanced into the aneurysm







Both balloons were then repositioned within the aneurysm, the right was inflated to protect the right PCA origin and also to keep the left balloon in the midline and slightly to the left within the aneurysm neck to protect the SCAs. It was decided that patient may tolerate left PCA sacrifice.







Multiple coils were placed through the simultaneously inflated Scepter balloons using the dual microcatheter technique.





Skull radiographs to show baseline coil mass and assess for future coil compaction.






76 year old with SAH




















What about What about recurrent aneurysms??


- Fill it with Pebbles not rocks
- Consider balloon or stent coiling or flow diverter
- No one said that you cannot Web it or Onyx it?!?
- I strongly prefer LVIS Blue

Durability of Rx: 63 year old patient with residual from previously coiled ruptured wide necked pericallosal aneurysm.




63 year old patient with residual from previously coiled ruptured wide necked pericallosal aneurysm.







63 year old patient with residual from previously coiled ruptured wide necked pericallosal aneurysm




3 x 15 Neuroform stent was placed and Duo catheter was used to access neck of aneurysm







Hydrocoils were used






Pre and post: 10 years on







40 year old with headaches and previously clipped aneurysm from the fronto-polar artery




Superselective WADA through Scepter




Balloon test occlusion with Scepter XC




Coiling through Scepter XC




Onyx through Scepter XC



Aneurysm no longer fills: Patient sent home the next day



Unruptured Aneurysms





Stent coiling with LVIS Blue: 45 year old with unruptured PCOM aneurysm







Stent assisted coiling after jailing and partial deployment





Compaction of stent tines at neck of aneurysm







Final DSA




 Significant ability to manipulate coverage of neck.


Flow diversion






Flow Diversion







PED + hydrocoils: 60 year old with right eye retro-ocular pain




- Giant partly thrombosed ICA communicating aneurysm.
- Decision to partly coil and placed PED.


DSA images from 3D






Placement of coils through Scepter XC






PED placement x 1






S/P PED x 2 and 18 months





Stent coiling with the LVIS Jr: 50 year old with growing ACOM aneurysm




Initial DSA





Road Map







Headway Duo 156





2.5 mm LVIS Jr in place







DSA post stent placement





Delayed Coil Embolization





Final result






Stent coiling with LVIS Jr: 58 year old with incidentally discovered 7 mm ACOM aneurysm





58 year old with incidentally discovered 7 mm ACOM aneurysm




DSA with 4 x 11 mm Scepter XC and Headway 17 in place






Balloon Assisted coiling






LVIS Jr placed through Scepter XC





Stent coiling with LVIS Jr: 56 year old with complex top of basilar SCA aneurysm






Placement of LVIS Jr in LEFT PCA






DSA S/P stent placement in left PCA






Y Stenting with second stent in Right PCA followed by coiling






Y Stenting with second stent in Right PCA followed by coiling





Final DSA result






LVIS Jr: Y stent: 35 year old with HIV






Stent coiling through the LVIS Jr






Pre and post







WEB + Coils: 63 year old with ruptured top of basilar aneurysm























Sofia in PCOM origin from ICA








43 year old with recurrence of previously clipped aneurysm














Fusiform Aneurysms: 75 y/o with Right P1 aneurysm





Left vertebral artery






5 F Sofia







LVIS Jr





Stent assisted coil embolization







Stent coil with LVIS Blue: 63 year old with ataxia and dysphagia





Giant partially thrombosed VBJ aneurysm






DSA with bilateral vertebral artery injection







Six overlapping 5.5 x 33 mm LVIS Blue Stents






Six overlapping 5.5 x 33 mm LVIS Blue Stents






Coil Embolization through Scepter XC







Final Construct








Post Embolization DSA








































































WEB treatment for distal intracranial aneurysms




82 year old with pericallosal aneurysm



















77 year old with PICA aneurysm






















WEB device for side wall aneurysms








4.1 mm (Geometric)


4.5 mm (Geometric)





















Salvage for intra-procedural rupture
































• Easy does it












































4 > H 00:00 01:17

ear old with lental ırysms

Thank you